


Essential Life Solutions has identified a need in laboratories to store columns efficiently when fully packed in a vertical condition. This is critical to the position of the column and the packed resin to remain vertical.

## Column Cart

We now proudly offer our COLUMN CART in three sizes to fit in any Lab setting depending on available space. Cart is constructed of chemically resistant aluminum and includes stainless steel casters with locking brakes. Columns can be mounted on both sides of the cart. Brackets can accommodate any columns from 10mm to 50mm diameter (any brand).



SIZE	ELS PART NUMBER	DIMENSIONS (in.) (L x W x H)	WEIGHT (lb.) EMPTY (approx.)
SMALL	ELS-COLUMN CART-SMALL	48 x 34 x 69	135
MEDIUM	ELS-COLUMN CART-MEDIUM	63 x 34 x 69	165
LARGE	ELS-COLUMN CART-LARGE	84 x 34 x 69	225



Column brackets (Blue) can be ordered separately to accommodate the appropriate amount of columns to be mounted.

*Column bracket assemblies sold separately,  
ELS PN: ELS-COLUMNBRACKETASSEMBLY*

**essentialLife Solutions Ltd**

308 Tosca Drive, Stoughton, MA 02072 USA

Phone 781-341-7240

Fax 781-341-7241

sales@essential-life.net

www.essential-life.net